

Recombinant ASFV g5Rp/D250R Protein, N-His

Catalog No: #AP97949

Package Size: #AP97949-1 50ug #AP97949-2 100ug #AP97949-3 1mg

Description

Product Name	Recombinant ASFV g5Rp/D250R Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	mRNA-decapping protein g5R, g5Rp, 3.1.3.-, ASFV-DP, Diphosphoinositol polyphosphate phosphohydrolase, DIPPP, 3.6.1.52, Ba71V-102, D250R
Accession No.	P32092
Target Species	African swine fever virus (ASFV)
Calculated MW	32.04 kDa
Target Length	Met1-His250
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.