

Recombinant ASFV Q706L Protein, N-His

Catalog No: #AP97951

Package Size: #AP97951-1 50ug #AP97951-2 100ug #AP97951-3 1mg

Description

Product Name	Recombinant ASFV Q706L Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Termination factor NPH-I homolog, 3.6.4.-, Ba71V-122, Q706L
Accession No.	Q89581
Target Species	African swine fever virus (ASFV)
Calculated MW	24.33 kDa
Target Length	Ile73-Gly267
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.