

Recombinant ASFV CP2475L (p37) Protein, N-His

Catalog No: #AP97952

Package Size: #AP97952-1 50ug #AP97952-2 100ug #AP97952-3 1mg

Description

Product Name	Recombinant ASFV CP2475L (p37) Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Polyprotein pp220, 220 kDa polyprotein, p34, p14, p37, p150, p5, Ba71V-92, CP2475L
Accession No.	Q08358
Target Species	African swine fever virus (ASFV)
Calculated MW	43.78 kDa
Target Length	Ala523-Gly893
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.