

Recombinant ASFV E248R/Ba71V-132 Protein, N-His

Catalog No: #AP97954

Package Size: #AP97954-1 50ug #AP97954-2 100ug #AP97954-3 1mg

Description

Product Name	Recombinant ASFV E248R/Ba71V-132 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Inner membrane protein pE248R, pE248R, Ba71V-132, E248R
Accession No.	Q65200
Target Species	African swine fever virus (ASFV)
Calculated MW	24.14 kDa
Target Length	Met1-Asn199
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.