## **Product Datasheet**

## Recombinant ASFV MGF\_360-4L Protein, N-His-SUMO

Catalog No: #AP97957

Package Size: #AP97957-1 50ug #AP97957-2 100ug #AP97957-3 1mg



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description	

Product Name	Recombinant ASFV MGF_360-4L Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MGF_360-4L
Accession No.	A0A385KMC2
Target Species	African swine fever virus (ASFV)
Calculated MW	25.02 kDa
Target Length	Asn14-Pro126
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store
	at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.