

Recombinant ASFV A224L Protein, N-His

Catalog No: #AP97959

Package Size: #AP97959-1 50ug #AP97959-2 100ug #AP97959-3 1mg

Description

Product Name	Recombinant ASFV A224L Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	A224L, A224L protein, Inhibitor of apoptosis protein, Protein A224L
Accession No.	A0A2X0THP1
Target Species	African swine fever virus (ASFV)
Calculated MW	28.80 kDa
Target Length	Met1-Leu224
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.