

Recombinant Pig CD163 Protein, C-His

Catalog No: #AP97965

Package Size: #AP97965-1 50ug #AP97965-2 100ug #AP97965-3 1mg

Description

Product Name	Recombinant Pig CD163 Protein, C-His
Host Species	Yeast
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Scavenger receptor cysteine-rich type 1 protein M130, CD163, Soluble CD163, sCD163, CD163, M130
Accession No.	Q2VL90
Target Species	Sus scrofa (Pig)
Calculated MW	14.76 kDa
Target Length	Arg232-Ser331
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.