

Recombinant CSFV E2 glycoprotein Protein, C-Strep-His

Catalog No: #AP97966

Package Size: #AP97966-1 50ug #AP97966-2 100ug #AP97966-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant CSFV E2 glycoprotein Protein, C-Strep-His
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E2 glycoprotein
Accession No.	B9VS89
Target Species	Classical swine fever virus (CSFV)
Calculated MW	42.19 kDa
Target Length	Arg1-Asp337
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.