

Recombinant Human CD221/IGF1R Protein, C-His

Catalog No: #AP97967

Package Size: #AP97967-1 50ug #AP97967-2 100ug #AP97967-3 1mg

Description

Product Name	Recombinant Human CD221/IGF1R Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	CD221, Insulin-like growth factor 1 receptor, IGF1R, Insulin-like growth factor I receptor, IGF-I receptor
Accession No.	P08069
Target Species	Homo sapiens (Human)
Calculated MW	56.37 kDa
Target Length	Met1-Val492
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.