

Recombinant Human MUC4 Protein, N-GST

Catalog No: #AP97969

Package Size: #AP97969-1 50ug #AP97969-2 100ug #AP97969-3 1mg

Description

Product Name	Recombinant Human MUC4 Protein, N-GST
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MUC-4, Ascites sialoglycoprotein 1, Ascites sialoglycoprotein 2, ASGP-1, ASGP-2, Testis mucin, MUC4, Ascites sialoglycoprotein, Pancreatic adenocarcinoma mucin, ASGP, Tracheobronchial mucin, Mucin-4
Accession No.	Q99102
Target Species	Homo sapiens (Human)
Calculated MW	30.79 kDa
Target Length	Gln5118-Phe5157
Tag Info	GST-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.