

Recombinant Streptococcus pyogenes M protein/emm12 Protein, N-His-KSI

Catalog No: #AP97973

Package Size: #AP97973-1 50ug #AP97973-2 100ug #AP97973-3 1mg

Description

Product Name	Recombinant Streptococcus pyogenes M protein/emm12 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	M protein, serotype 12, emm12
Accession No.	P19401
Target Species	Streptococcus pyogenes
Calculated MW	28.14 kDa
Target Length	Asp42-Leu91
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.