

Recombinant Streptococcus pneumoniae PBP3 Protein, C-His

Catalog No: #AP97979

Package Size: #AP97979-1 50ug #AP97979-2 100ug #AP97979-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Streptococcus pneumoniae PBP3 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	D-alanyl-D-alanine carboxypeptidase PBP3, pbp3
Accession No.	WP_061743357.1
Target Species	Streptococcus pneumoniae
Calculated MW	44.11 kDa
Target Length	Gln23-Leu413
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.