

Recombinant Streptococcus agalactiae CYK3/IgdE3 Protein, N-GST & C-His

Catalog No: #AP97983

Package Size: #AP97983-1 50ug #AP97983-2 100ug #AP97983-3 1mg

Description

Product Name	Recombinant Streptococcus agalactiae CYK3/IgdE3 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	hypothetical protein SAG0135_08530, hypothetical protein SAG0135_08530
Accession No.	EPU23939.1
Target Species	Streptococcus agalactiae
Calculated MW	76.72 kDa
Target Length	Thr45-Tyr463
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.