

Recombinant Streptococcus agalactiae GBS Protein, C-His

Catalog No: #AP97988

Package Size: #AP97988-1 50ug #AP97988-2 100ug #AP97988-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Streptococcus agalactiae GBS Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	alpha-like surface protein, partial, alpha-like surface protein, partial
Accession No.	WP_133256377.1 & WP_148547943.1
Target Species	Streptococcus agalactiae
Calculated MW	39.73 kDa
Target Length	[Ala56-Pro229] & [Ser57-Ser226]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.