

Recombinant Streptococcus pyogenes EndoS Protein, N-GST

Catalog No: #AP97993

Package Size: #AP97993-1 50ug #AP97993-2 100ug #AP97993-3 1mg

Description

Product Name	Recombinant Streptococcus pyogenes EndoS Protein, N-GST
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Endo-beta-N-acetylglucosaminidase
Accession No.	Q9APG4
Target Species	Streptococcus pyogenes
Calculated MW	134.61 kDa
Target Length	Glu37-Lys995
Tag Info	GST-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.