

Recombinant JEV NS3/Serine protease NS3 Protein, N-His

Catalog No: #AP98003

Package Size: #AP98003-1 50ug #AP98003-2 100ug #AP98003-3 1mg

Description

Product Name	Recombinant JEV NS3/Serine protease NS3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Serine protease NS3, 3.4.21.91, 3.6.1.15, 3.6.4.13, Flavivirin protease NS3 catalytic subunit, Non-structural protein 3
Accession No.	P27395
Target Species	Japanese encephalitis virus (strain SA-14) (JEV)
Calculated MW	51.69 kDa
Target Length	Pro1685-Arg2123
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.