

# Recombinant *Neisseria meningitidis* pilE Protein, N-His

Catalog No: #AP98006

Package Size: #AP98006-1 50ug #AP98006-2 100ug #AP98006-3 1mg

## Description

|                |                                                                                                                                                                           |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant <i>Neisseria meningitidis</i> pilE Protein, N-His                                                                                                             |
| Host Species   | <i>E. coli</i>                                                                                                                                                            |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                           |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| Other Names    | Pilin, pilE                                                                                                                                                               |
| Accession No.  | Q51237                                                                                                                                                                    |
| Target Species | <i>Neisseria meningitidis</i>                                                                                                                                             |
| Calculated MW  | 17.10 kDa                                                                                                                                                                 |
| Target Length  | Gln25-Ser163                                                                                                                                                              |
| Tag Info       | His-tag                                                                                                                                                                   |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.