

Recombinant *Neisseria meningitidis* neuA Protein, C-His-Strep

Catalog No: #AP98014

Package Size: #AP98014-1 50ug #AP98014-2 100ug #AP98014-3 1mg

Description

Product Name	Recombinant <i>Neisseria meningitidis</i> neuA Protein, C-His-Strep
Host Species	<i>E. coli</i>
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	N-acetylneuraminate cytidyltransferase, 2.7.7.43, CMP-N-acetylneuraminic acid synthase, CMP-NeuNAc synthase, CMP-sialic acid synthase, neuA, siaB, synB
Accession No.	P0A0Z8
Target Species	<i>Neisseria meningitidis</i>
Calculated MW	26.90 kDa
Target Length	Met1-Ser228
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.