

Recombinant Mouse S100A14 Protein, N-GST & C-His

Catalog No: #AP98015

Package Size: #AP98015-1 50ug #AP98015-2 100ug #AP98015-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse S100A14 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	S100A14, S100 calcium-binding protein A14, S100A15, S114, Protein S100-A14
Accession No.	Q9D2Q8
Target Species	Mus musculus (Mouse)
Calculated MW	38.51 kDa
Target Length	Met1-Lys93
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.