

Recombinant MPXV H3L Protein, C-His

Catalog No: #AP98018

Package Size: #AP98018-1 50ug #AP98018-2 100ug #AP98018-3 1mg

Description

Product Name	Recombinant MPXV H3L Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	H3L, H3L
Accession No.	NP_536520.1
Target Species	Monkeypox virus (strain Zaire-96-I-16) (MPX)
Calculated MW	33.92 kDa
Target Length	Ala4-Gly270
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.