

Recombinant HAdV-41 E1B-55K Protein, N-His

Catalog No: #AP98021

Package Size: #AP98021-1 50ug #AP98021-2 100ug #AP98021-3 1mg

Description

Product Name	Recombinant HAdV-41 E1B-55K Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E1B 55 kDa protein, E1B-55K, E1B protein, large T-antigen, E1B-495R
Accession No.	P10546
Target Species	Human adenovirus F serotype 41 (HAdV-41) (Human adenovirus 41)
Calculated MW	44.99 kDa
Target Length	Lys98-Asp472
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.