

Recombinant HAdV-5 Protease/L3 Protein, C-His

Catalog No: #AP98022

Package Size: #AP98022-1 50ug #AP98022-2 100ug #AP98022-3 1mg

Description

Product Name	Recombinant HAdV-5 Protease/L3 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Protease, 3.4.22.39, Adenain, Adenovirus protease, AVP, Adenovirus proteinase, Endoprotease, L3
Accession No.	P03253
Target Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Calculated MW	24.13 kDa
Target Length	Met1-Met204
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.