

Recombinant HAdV-3 Protease/L3 Protein, C-His

Catalog No: #AP98024

Package Size: #AP98024-1 50ug #AP98024-2 100ug #AP98024-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HAdV-3 Protease/L3 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Protease, 3.4.22.39, Adenain, Adenovirus protease, AVP, Adenovirus proteinase, Endoprotease, L3
Accession No.	P10381
Target Species	Human adenovirus B serotype 3 (HAdV-3) (Human adenovirus 3)
Calculated MW	24.86 kDa
Target Length	Met1-Gln209
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.