

Recombinant HAdV-5 Penton protein/L2 Protein, C-Strep

Catalog No: #AP98025

Package Size: #AP98025-1 50ug #AP98025-2 100ug #AP98025-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HAdV-5 Penton protein/L2 Protein, C-Strep
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Penton protein, CP-P, Penton base protein, Protein III, L2
Accession No.	P12538
Target Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Calculated MW	59.48 kDa
Target Length	Thr49-Phe571
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.