

Recombinant HAdV-5 Fiber protein/L5 Protein, N-His

Catalog No: #AP98028

Package Size: #AP98028-1 50ug #AP98028-2 100ug #AP98028-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HAdV-5 Fiber protein/L5 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Fiber protein, SPIKE, Protein IV, L5
Accession No.	P11818
Target Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Calculated MW	22.52 kDa
Target Length	Gly387-Glu581
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.