

Recombinant HBV-B HBsAg Protein, N-His

Catalog No: #AP98029

Package Size: #AP98029-1 50ug #AP98029-2 100ug #AP98029-3 1mg

Description

Product Name	Recombinant HBV-B HBsAg Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Large envelope protein, L glycoprotein, L-HBsAg, LHB, Large S protein, Large surface protein, Major surface antigen, S
Accession No.	Q67926
Target Species	Hepatitis B virus genotype B2 subtype adw (isolate China/patient4/1996) (HBV-B)
Calculated MW	29.11 kDa
Target Length	Met1-Arg253
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.