

Recombinant Human Immunoglobulin G3 (IgG3)

Catalog No: #AP98031

Package Size: #AP98031-1 50ug #AP98031-2 100ug #AP98031-3 1mg

Description

Product Name	Recombinant Human Immunoglobulin G3 (IgG3)
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Immunoglobulin G3, IgG3
Accession No.	P01860
Target Species	Homo sapiens (Human)
Calculated MW	150 kDa
Target Length	Full length.
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.