

Recombinant Aspergillus flavus Uricase/Uox Protein, C-hFc (IgG1)

Catalog No: #AP98056

Package Size: #AP98056-1 50ug #AP98056-2 100ug #AP98056-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Aspergillus flavus Uricase/Uox Protein, C-hFc (IgG1)
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Uricase, 1.7.3.3, Rasburicase, Urate oxidase, uaZ, Uox
Accession No.	Q00511
Target Species	Aspergillus flavus
Calculated MW	62.45 kDa
Target Length	Met1-Lys302
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.