

Recombinant Human SMC1A Protein, N-His

Catalog No: #AP98058

Package Size: #AP98058-1 50ug #AP98058-2 100ug #AP98058-3 1mg

Description

Product Name	Recombinant Human SMC1A Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	DXS423E, SB1.8, KIAA0178, Structural maintenance of chromosomes protein 1A, SMC-1A, Sb1.8, SMC1, SMC1A, SMC protein 1A, SMC1L1, SMC-1-alpha
Accession No.	Q14683
Target Species	Homo sapiens (Human)
Calculated MW	22.58 kDa
Target Length	Ile489-Asp667
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.