

Recombinant PEDV S/Spike glycoprotein Protein, C-Strep

Catalog No: #AP98061

Package Size: #AP98061-1 50ug #AP98061-2 100ug #AP98061-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant PEDV S/Spike glycoprotein Protein, C-Strep
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, S
Accession No.	A0A059VFK8
Target Species	Porcine epidemic diarrhea virus
Calculated MW	49.65 kDa
Target Length	Pro20-Ile461
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.