

Recombinant Human Osteocalcin/BGLAP/OC Protein, N-GST & C-His

Catalog No: #AP98064

Package Size: #AP98064-1 50ug #AP98064-2 100ug #AP98064-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Osteocalcin/BGLAP/OC Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Gamma-carboxyglutamic acid-containing protein, BGLAP, BGP, Bone Gla protein, Osteocalcin
Accession No.	P02818
Target Species	Homo sapiens (Human)
Calculated MW	36.94 kDa
Target Length	Lys24-Val100
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.