

Recombinant Mouse CD74 Protein, C-His

Catalog No: #AP98066

Package Size: #AP98066-1 50ug #AP98066-2 100ug #AP98066-3 1mg

Description

Product Name	Recombinant Mouse CD74 Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	H-2 class II histocompatibility antigen gamma chain, Ia antigen-associated invariant chain, Ii, Ii chain, MHC class II-associated invariant chain, CD74, Class-II-associated invariant chain peptide, CLIP, Cd74
Accession No.	P04441
Target Species	Mus musculus (Mouse)
Calculated MW	28.34 kDa
Target Length	Gln56-Leu279
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.