

Recombinant Human METTL5 Protein, N-His

Catalog No: #AP98075

Package Size: #AP98075-1 50ug #AP98075-2 100ug #AP98075-3 1mg

Description

Product Name	Recombinant Human METTL5 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Methyltransferase-like protein 5, METTL5, rRNA N6-adenosine-methyltransferase METTL5
Accession No.	Q9NRN9
Target Species	Homo sapiens (Human)
Calculated MW	25.88 kDa
Target Length	Met1-Phe209
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.