

Recombinant Mouse CD257/BAFF/TNFSF13B Protein, N-Fc

Catalog No: #AP98081

Package Size: #AP98081-1 50ug #AP98081-2 100ug #AP98081-3 1mg

Description

Product Name	Recombinant Mouse CD257/BAFF/TNFSF13B Protein, N-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tumor necrosis factor ligand superfamily member 13B, B-cell-activating factor, BAFF, CD257, Tumor necrosis factor ligand superfamily member 13b, membrane form, Tumor necrosis factor ligand superfamily member 13b, soluble form, Tnfsf13b, Baff
Accession No.	Q9WU72
Target Species	Mus musculus (Mouse)
Calculated MW	48.44 kDa
Target Length	Arg126-Leu309
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.