

Recombinant SFTSV NP/Nucleoprotein Protein, N-His

Catalog No: #AP98087

Package Size: #AP98087-1 50ug #AP98087-2 100ug #AP98087-3 1mg

Description

Product Name	Recombinant SFTSV NP/Nucleoprotein Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	nucleocapsid protein, nucleocapsid protein
Accession No.	YP_006504092.1
Target Species	SFTS phlebovirus (isolate SFTSV/Human/China/HB29/2010) (Severe fever with thrombocytopenia virus)
Calculated MW	31.41 kDa
Target Length	Met1-Leu245
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.