

Recombinant MARV GP2/Envelope glycoprotein 2 Protein, C-Fc

Catalog No: #AP98095

Package Size: #AP98095-1 50ug #AP98095-2 100ug #AP98095-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant MARV GP2/Envelope glycoprotein 2 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Envelope glycoprotein, GP1,2, GP, Virion spike glycoprotein, GP2
Accession No.	Q1PDC7
Target Species	Lake Victoria marburgvirus (strain Ravn-87) (MARV) (Marburg virus (strain Kenya/Ravn/1987))
Calculated MW	52.61 kDa
Target Length	Ser436-Trp649
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.