

Recombinant CSFV NS5B Protein, C-His

Catalog No: #AP98098

Package Size: #AP98098-1 50ug #AP98098-2 100ug #AP98098-3 1mg

Description

Product Name	Recombinant CSFV NS5B Protein, C-His
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	polyprotein, polyprotein, NS5B, RNA-directed RNA polymerase, p75
Accession No.	AAT85717
Target Species	Classical swine fever virus (strain Brescia) (CSFV) (Hog cholera virus)
Calculated MW	83.25 kDa
Target Length	Ser3181-Gly3897
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.