

Recombinant CSFV Glycoprotein E2/gp55 Protein, N-His

Catalog No: #AP98099

Package Size: #AP98099-1 50ug #AP98099-2 100ug #AP98099-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant CSFV Glycoprotein E2/gp55 Protein, N-His
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Envelope glycoprotein E2, gp55
Accession No.	P21530
Target Species	Classical swine fever virus (strain Brescia) (CSFV) (Hog cholera virus)
Calculated MW	41.62 kDa
Target Length	Ala739-Ala1029
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.