

Recombinant BVDV Glycoprotein E1/gp33 Protein, N-His

Catalog No: #AP98103

Package Size: #AP98103-1 50ug #AP98103-2 100ug #AP98103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant BVDV Glycoprotein E1/gp33 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Envelope glycoprotein E1, gp33
Accession No.	P19711
Target Species	Bovine viral diarrhea virus (isolate NADL) (BVDV) (Mucosal disease virus)
Calculated MW	20.97 kDa
Target Length	Ala498-Thr663
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.