

Recombinant BVDV Glycoprotein E2/gp55 Protein, C-Fc

Catalog No: #AP98104

Package Size: #AP98104-1 50ug #AP98104-2 100ug #AP98104-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant BVDV Glycoprotein E2/gp55 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E2 glycoprotein, E2 glycoprotein, gp55
Accession No.	AGC12054.1
Target Species	Bovine viral diarrhea virus (BVDV) (Mucosal disease virus)
Calculated MW	52.78 kDa
Target Length	Cys1-Asp218
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.