

Recombinant BVDV Glycoprotein E2/gp53 Protein, C-His

Catalog No: #AP98105

Package Size: #AP98105-1 50ug #AP98105-2 100ug #AP98105-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant BVDV Glycoprotein E2/gp53 Protein, C-His
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E2 glycoprotein, E2 glycoprotein, gp53
Accession No.	Q65808
Target Species	Bovine viral diarrhea virus (BVDV) (Mucosal disease virus)
Calculated MW	43.07 kDa
Target Length	His36-His370
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.