

Recombinant EEEV Spike glycoprotein E1 Protein, N-His

Catalog No: #AP98107

Package Size: #AP98107-1 50ug #AP98107-2 100ug #AP98107-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant EEEV Spike glycoprotein E1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Structural polyprotein, p130, Spike glycoprotein E1, E1 envelope glycoprotein
Accession No.	P08768
Target Species	Eastern equine encephalitis virus (EEEV) (Eastern equine encephalomyelitis virus)
Calculated MW	48.07 kDa
Target Length	Tyr799-Gly1215
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.