

Recombinant MACV GP1 Protein, C-Fc

Catalog No: #AP98112

Package Size: #AP98112-1 50ug #AP98112-2 100ug #AP98112-3 1mg

Description

Product Name	Recombinant MACV GP1 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	glycoprotein precursor, glycoprotein precursor, GP1, Glycoprotein G1
Accession No.	NP_899212.1
Target Species	Machupo virus (MACV)
Calculated MW	48.56 kDa
Target Length	Glu87-Leu261
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.