

## Recombinant Human MNAT1 Protein, N-His

Catalog No: #AP98116

Package Size: #AP98116-1 50ug #AP98116-2 100ug #AP98116-3 1mg

## Description

|                |                                                                                                                                                                                                                    |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Human MNAT1 Protein, N-His                                                                                                                                                                             |
| Host Species   | E. coli                                                                                                                                                                                                            |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                                                                    |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                                                                    |
| Other Names    | p36, MNAT1, MAT1, Menage a trois, Cyclin-G1-interacting protein, RING finger protein MAT1, RING finger protein 66, RNF66, CDK7/cyclin-H assembly factor, p35, CAP35, CDK-activating kinase assembly factor<br>MAT1 |
| Accession No.  | P51948                                                                                                                                                                                                             |
| Target Species | Homo sapiens (Human)                                                                                                                                                                                               |
| Calculated MW  | 14.62 kDa                                                                                                                                                                                                          |
| Target Length  | Met1-Ile113                                                                                                                                                                                                        |
| Tag Info       | His-tag                                                                                                                                                                                                            |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                                         |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.                                              |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.