

Recombinant Human KANK2 Protein, N-His

Catalog No: #AP98118

Package Size: #AP98118-1 50ug #AP98118-2 100ug #AP98118-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KANK2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ankyrin repeat domain-containing protein 25, SRC-interacting protein, SRC-1-interacting protein, SRC1-interacting protein, ANKRD25, KIAA1518, MXRA3, KN motif and ankyrin repeat domain-containing protein 2, Matrix-remodeling-associated protein 3, KANK2, SIP
Accession No.	Q63ZY3
Target Species	Homo sapiens (Human)
Calculated MW	21.00 kDa
Target Length	Met1-Asn161
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.