

Recombinant Human RBBP9 Protein, N-His

Catalog No: #AP98120

Package Size: #AP98120-1 50ug #AP98120-2 100ug #AP98120-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RBBP9 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	BOG, RBBP9, B5T-overexpressed gene protein, RBBP10, Serine hydrolase RBBP9, RBBP-9, Retinoblastoma-binding protein 9, Protein BOG, Retinoblastoma-binding protein 10, RBBP-10
Accession No.	O75884
Target Species	Homo sapiens (Human)
Calculated MW	20.63 kDa
Target Length	Gly66-Glu231
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.