

Recombinant Mouse SIGLEC12 Protein, N-His

Catalog No: #AP98122

Package Size: #AP98122-1 50ug #AP98122-2 100ug #AP98122-3 1mg

Description

Product Name	Recombinant Mouse SIGLEC12 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Sialic acid-binding Ig-like lectin 12, Siglec-L1, Sialic acid-binding Ig-like lectin-like 1, SIGLEC12, SLG, Siglec-12, SIGLECL1
Accession No.	Q91Y57
Target Species	Mus musculus (Mouse)
Calculated MW	28.98 kDa
Target Length	Gln19-Ala253
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.