

Recombinant Human TNFRSF25 Protein, C-Fc

Catalog No: #AP98124

Package Size: #AP98124-1 50ug #AP98124-2 100ug #AP98124-3 1mg

Description

Product Name	Recombinant Human TNFRSF25 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Apoptosis-inducing receptor AIR, Apoptosis-mediating receptor DR3, Apoptosis-mediating receptor TRAMP, Protein WSL-1, LARD, Lymphocyte-associated receptor of death, Tumor necrosis factor receptor superfamily member 25, Apo-3, Protein WSL, APO3, DR3, TNFRSF12, WSL, TNFRSF25, Death receptor 3, DDR3, WSL1
Accession No.	Q93038
Target Species	Homo sapiens (Human)
Calculated MW	47.74 kDa
Target Length	Met1-Gln199
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.