

Recombinant Human PDE4B Protein, N-His

Catalog No: #AP98132

Package Size: #AP98132-1 50ug #AP98132-2 100ug #AP98132-3 1mg

Description

Product Name	Recombinant Human PDE4B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	PDE4B, cAMP-specific 3',5'-cyclic phosphodiesterase 4B, PDE32, DPDE4
Accession No.	Q07343
Target Species	Homo sapiens (Human)
Calculated MW	39.32 kDa
Target Length	Asn334-Ser654
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.