

Recombinant Human PI3 Protein, N-GST & C-His

Catalog No: #AP98133

Package Size: #AP98133-1 50ug #AP98133-2 100ug #AP98133-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PI3 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	PI3, Elastase-specific inhibitor, WAP3, Elafin, Protease inhibitor WAP3, Peptidase inhibitor 3, PI-3, Skin-derived antileukoprotease, SKALP, ESI, WFDC14, WAP four-disulfide core domain protein 14
Accession No.	P19957
Target Species	Homo sapiens (Human)
Calculated MW	38.12 kDa
Target Length	Ala23-Gln117
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.